

Clinispin® CT25 Centrifuge

with Microprocessor & Brushless DC Motor

The essential centrifuge for Blood Processing in
Clinics/Hospitals, Veterinary Surgeries &
Contract Research Organisations



**LIGHTWEIGHT, PORTABLE
& MAINTENANCE FREE**



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

Clinispin® CT25 Centrifuge

with Microprocessor & Brushless DC Motor

The Clinispin CT25 is a lightweight, portable, maintenance free centrifuge designed for blood centrifugation and other clinical applications, including PRP.

FEATURES

- Variable speed with precise RPM/RCF setting from 500 to 6500 RPM with least count of 100 RPM
- Lightweight and portable
- Universal voltage and maintenance free
- Accommodates most common collection tubes - comes with adaptors for 5ml, 7ml & 1.5/2ml tubes
- Supplied with 8 x 15ml rotor
- One switch change over from RPM to RCF
- Smart air flow design ensures quiet operation and low temperature rise
- Automatic imbalance detection for safe operation



SPECIFICATIONS

Fixed Angle Rotor: 8 x 15ml (Autoclavable) - 3873g

Max. RPM: 6500

Motor Type: Brushless DC

Speed Setting: 500-6500 RPM

Timer Setting: 1-30 minutes and infinite mode

Lid Lock Safety: Yes

Noise Level: <60dB

Imbalance Detection: Yes

Dimensions: 260 x 244 x 205 mm

Weight: 4kg

Adaptors (included in standard pack): Chosen rotor with (adaptors for 5ml, 7ml and 1.5/2ml tubes - takes most standard blood collection tubes)

SUGGESTED APPLICATIONS

- Clinics/Hospitals/General Practitioners
- Veterinary Surgeries
- Contract Research/Pharmaceutical Organisations
- Laboratories
- Remote/Rural Healthcare Settings



8 x 15ml (Autoclavable)



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**