

Clinispin™ Mini

Mini Veterinary Centrifuge



**SMALL FOOTPRINT,
FOR SMALL WORKPLACES**

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



Mini Veterinary Centrifuge

The Clinispin Mini offers versatility, simplicity, safety and ease of use in a single centrifuge. It is suitable for microfiltration, cell separation, centrifugal tube sedimentation and spin down applications.

FEATURES

- 8 position mini centrifuge tube rotor
- Digital display with microprocessor control
- Set variable speed and timer
- Electronic safety brake on lid opening for user safety
- Magnetic rotor for easy tool free changing of rotor
- Starts and stops in seconds
- Lid safety interlock
- Smooth and silent, discreet running
- Small footprint, can fit in small workplaces
- Interchangeable magnetic rotors (MCT rotor & PCR rotor)
- Low voltage power consuming device



SPECIFICATIONS

Motor Type: Brushless DC Motor

Rotor Capacity: 8 x 1.5/2.0ml Tube Rotor & 2 Place PCR Strip Rotor

Speed: 1000-6000 RPM, 2000 xg (Max. RCF)

Speed Accuracy: $\pm 2\%$ RPM

Run Time & Mode: 1-25 mins & infinite mode

Dimensions: 162 (W) x 157 (D) x 116 (H) mm

Weight: 1.1kg

Power Adapter: Input: 100-240 VAC, 50/60 Hz, Output: 24V \equiv 1.5A

Power Consumption: 15W

Permissible Relative Humidity: $\leq 80\%$

Ambient Temperature: +5°C to +40°C

Altitude: Use up to an altitude of 2000m above MSL (mean sea level)

Pollution Degree: 2

Environment: For indoor use only



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**