

www.woodleyequipment.com

g-Pet

GLUCOSE
Veterinary Monitoring System



**Making Testing
Quick, Easy &
Economical!**

Complete Kit Includes:

- Meter**
- Test strips**
- Control**
- User Guide**
- Log Book**

Tel: +44 (0)1204 669033
Fax: +44 (0)1204 669034
Email: sales@woodleyequipment.com
Web www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE. BL6 7NY.

Making Testing Easy!



Step 1

Insert the test strip
(Meter will turn on automatically)



Step 2

Apply blood



Step 3

Results in 10 seconds

g-Pet

GLUCOSE
Veterinary Monitoring System

- Biosensor auto code technology
- 1.5 microlitre blood sample
- 4 reminder alarms for worry free testing
- 8 second countdown
- 300 result memory with date and time
- Accurate results
- Big LCD screen
- Easy to use
- Hypoglycaemic and hyperglycaemic indicator alarms
- PC download option

Tel: +44 (0)1204 669033

Fax: +44 (0)1204 669034

Email: sales@woodleyequipment.com

Web www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE. BL6 7NY.