

www.woodleyequipment.com

nova Vet

STATSENSOR XPRESS LACTATE METER



- Results in 13 seconds
- Easy to use
- Portable
- Small 0.6µl sample size
- Capillary, venous or arterial whole blood
- Stores 400 Patient results
- Corrects result for HCT
- No user calibration
- Eliminates interfering substances

NovaVet StatSensor Xpress Lactate Meter for use in the quantitative determination of Lactate in fresh capillary, venous and arterial whole blood.

Hand-held, battery-powered, in vitro diagnostic laboratory instrument that works in conjunction with Nova Biomedical lactate electrochemical test strips to measure lactate in a whole blood sample.

The meter can store up to 400 patient and/or quality control test results and has a built-in beeper for audible alerts and prompts.

Tel: +44 (0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

WOODLEY
EQUIPMENT COMPANY LTD.

www.woodleyequipment.com

nova
biomedical



*NovaVet
High Impact Carry Case*

(Case and contents for illustration only)

Specification:

Temperature: 59°F-104°F (15°C- 40°C)

Altitude: Up to 15,000 feet

Humidity: 10% to 90% relative humidity

Weight: 75g (2.65 oz)

Size: 91.4 x 58.4 x 22.9 mm

Lactate Measurement Range:

0.3-20.0 mmol/L

2.7-180 mg/dL

Sample Types:

Whole Blood: Capillary, Arterial, Venous

Tel: +44 (0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

WOODLEY
EQUIPMENT COMPANY LTD.