

www.woodleyequipment.com

ANNOUNCING

THE NEXT GENERATION OF HANDHELD CRITICAL CARE ANALYSERS



epoc
BLOOD ANALYSIS

:: a cost effective, state of the art point of care solution ::

The epoc, developed by the founder of the i-STAT analyser, is the first wireless point of care testing solution to use Smart Card technology. The epoc provides immediate blood gas, electrolyte and critical care biochemistry results directly to a handheld PDA (included) and wirelessly interfaces to an electronic data manager (EDM) and bluetooth printer. It is a completely affordable analyser that is easy to use and requires no refrigeration of consumables.

**NOW YOU CAN UPGRADE YOUR
PORTABLE CRITICAL CARE ANALYSER
TO THE VERY LATEST TECHNOLOGY**

Contact us for a demonstration
and to discuss our **not to be
missed trade-in offer** for your
i-STAT® system*

* i-STAT® is a registered trade mark of Abbott Point of Care Inc

Tel: +44 (0)1204 669033
Fax: +44 (0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS . ST JOHNS STREET . HORWICH . BOLTON . LANCASHIRE . BL6 7NY

Key Features:

- Accurate blood gas, electrolyte and critical care biochemistry results in 30 seconds
- Room temperature cartridge storage
- 100 µl sample volume, eliminating overfill issues
- Wireless system that interfaces to an electronic data manager (EDM) and Bluetooth printer
- Veterinary specific software – does not require biannual software updates
- Easy to use with touch screen technology
- Completely portable handheld system with rechargeable battery
- Ready to use anytime, anywhere – in-practice, theatre, client site or yard
- Integrated calibration and quality checks
- Bar-coded smart cards for error free identification
- 2000 patient results stored on analyser
- Ability to input all patient data, temperature and even ventilator settings
- Cartridge batch numbers automatically linked to patient results
- New parameters continually being developed

Lower System and Consumable Costs than other Leading Portable Critical Care Analysers including the i-STAT® system*



epoc
BLOOD ANALYSIS

3 Simple Steps TO ACCURATE & IMMEDIATE RESULTS



:: state of the art blood analysis ::

Exclusively from Woodley Equipment Company Ltd

* i-STAT® is a registered trade mark of Abbott Point of Care Inc

Tel: +44 (0)1204 669033
Fax: +44 (0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS . ST JOHNS STREET . HORWICH . BOLTON . LANCASHIRE . BL6 7NY