

Digital Incubator 10 Litre

Available for hire or purchase for your clinical trial.



Overview

This 10L digital incubator is a compact design and is an economical and space saving solution for microbiology or haematology applications. The temperature can be freely adjusted in increments of 0.1°C up to a maximum of 70°C. Heating elements are incorporated into the housing base and walls to ensure temperature stability and uniformity.

The transparent window allows full visibility of samples without opening the door, thus preventing heat loss. The housing and door frame are constructed entirely of epoxy coated mild steel.

Operating Features

- Compact design
- Power light illuminated when switched on
- LED temperature display
- Excellent temperature stability
- Up and down key for temperature adjustment
- Supplied with one perforated stainless steel shelf
- Easy to use, safe and reliable

Specification

- Net Capacity: 10 Litres
- Temperature Range: Ambient +5°C to +70°C
- Temperature Accuracy: at 37°C +/- 0.5°C
- Temperature Stability: at 37°C +/- 0.5°C
- Controls: PID Control
- Power Supply: 230V, 50/60Hz, 0.3A
- Weight: 8.3kgs
- Internal Dimensions: 200h x 230w x 220d mm
- External Dimensions: 335h x 285w x 280d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Digital Incubator 10 Litre

Product Code: WD7INC001

Version: 10/19 V2

Product specifications can be subject to change without prior notice

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company