

# MiniSpin Plus Microcentrifuge

Available for hire or purchase for your clinical trial.



## Overview

The MiniSpin Plus personal centrifuge has a sleek exterior and features reliable and durable technology. The design combines functionality, safety, ease of use and comfort.

The MiniSpin Plus has 12-positions for research with small sample requirements. The perfectly designed metal rotor lid even accommodates virtually all spin columns.

The maintenance-free drive of the MiniSpin Plus accelerates up to 12 x 1.5/2.0 ml microcentrifuge tubes to a maximum speed of 14,100 x g (14,500 RPM).

## Operating Feature

- Exceptionally small footprint
- Extremely quiet - low noise levels
- Maintenance-free drive
- Easy-to-follow digital display
- Low sample heating (only 12°C / 53.6°F after 20 min. at max. speed)
- Separate short-spin key
- Automatic lid release
- Built-in power unit
- Autoclavable rotor (121°C / 250°F, 20 min.)
- Acceleration and braking times <13 sec. (at max. speed)
- Stainless steel rotor lid
- Includes refined "Black Line" aluminum rotor
- Max. speed of 14,500 RPM; max. RCF 14,100 x g
- RPM/RCF setting as required
- Continuous centrifugation possible
- Optional: adaptor for 0.2 ml PCR tubes, adaptors for 0.4 microcentrifuge tubes, 0.5 ml microcentrifuge tubes and 0.6 ml Microtainer®

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Rotor: 12 x 1.5 ml / 2.0 ml
- Max. RCF: 14,100 x g
- Max. Speed: 14,500 RPM
- Acceleration Time: 13 seconds
- Braking Time: 12 seconds
- Electronic Timer: 15 seconds to 99 minutes with continuous function
- Weight: 4.3 kg / 9.5 lbs
- Dimensions: 119h x 226w x 239d mm  
4.7h x 8.9w x 9.4d inches
- Voltage: 110V / 230V
- Frequency: 50/60 Hz



**Please contact us for a quotation for the hire or purchase of this product**

Product description: MiniSpin Plus microcentrifuge

Product Code: WD7CEN039

Version: V3 10/19

Product specifications can be subject to change without prior notice

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company