

Thermo Scientific™ Indiko™ Plus Veterinary Chemistry Analyser

With optional ISE module for rapid, accurate electrolyte analysis

A comprehensive system solution for veterinary testing offering maximum ease of use, reliability and cost-effectiveness for clinical chemistry and electrolyte assays.



Woodley Equipment Company Ltd represents
Thermo Scientific™ Indiko™ vet solution in the UK:
Trusted Partners in Veterinary Diagnostics



thermo
scientific

Authorised Distributor

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

Thermo Scientific™ Indiko™ Plus Veterinary Chemistry Analyser

Woodley Equipment Company Ltd represents Thermo Scientific™ Indiko™
vet solution in the UK: **Trusted Partners in Veterinary Diagnostics**

The Indiko Plus Veterinary Chemistry Analyser from Thermo Scientific is ideal for multi-tasking vet clinics and labs with a wide range of test requests. It can run multiple tests concurrently, with up to 108 samples and 42 reagents.

FEATURES

- Fully automated, random-access benchtop analyser for labs of any size
- Complete system solution: easy to use analyser, reagents, calibrators & consumables
- Accurate, reliable results within minutes
- Small sample volume (2-120µl)
- Can run multiple tests concurrently, with up to 108 samples and 42 reagents
- Optional integrated ISE module for rapid, accurate electrolyte analysis
- STAT mode for urgent samples
- Automated testing features reduces risk of errors
- Continuous access to samples, reagents and cuvettes ensures uninterrupted test runs
- Bidirectional LIS connection
- Low maintenance



thermo
scientific

Authorised Distributor

ASSAYS

- Alkaline Phosphatase (IFCC) Plus
 - ALT/GPT (IFCC)
 - Amylase (IFCC)
 - AST/GOT (IFCC)
 - Bicarbonate
 - Bilirubin direct
 - Bilirubin total (DCA)
 - Calcium
 - Cholesterol
 - Cholinesterase
 - CK (IFCC)
 - Creatinine (Enzymatic)
 - Creatinine (Jaffe)
 - Gamma-GT (IFCC)
 - Glucose (HK)
 - HDL-Cholesterol
 - Iron
 - LDH (IFCC)
 - LDL-Cholesterol
 - Lipase
 - Magnesium
 - Phosphorus
 - Pyridoxal Phosphate
 - Total Protein Plus
 - Triglycerides
 - U/CSF Protein
 - Urea
 - Uric Acid (AOX)
 - Glucose (GOD-POD)
- Open channels for speciality tests**

SPECIFICATIONS

Capacity: Up to 350 Photometric Tests/Hour
Up to 135 ISE Tests/Hour

On-board Sample Capacity: Max. 108 with 6 x 18-pos. sample rack

On-board Reagent Capacity: Max. 42 reagent positions

Sample Volumes: 2-120µl

Reagent Volumes: 2-240µl

Walk-away Time: Up to 3 hours

Water Consumption: 2.5 litres/hour

LIS Interface: CLSI LIS02-A2

Power Requirements: 100-240 VAC ±10%, 50-60 Hz ±5%, 300 W

Operating Temperature: +18°C to +30°C

Operating Humidity: 40% to 80% (non-condensing)

Dimensions (W x D x H): 94 x 70 x 62 cm/130cm with the open cover

Weight: 110kg



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**