

InSight[®] Canine-IA

Canine Immunoassay Analyser

A portable, handheld, easy to use Canine Immunoassay Analyser providing quick, accurate and reliable results.



TESTS AVAILABLE

- cPL
- T4
- TSH
- Cortisol
- Progesterone
- Alpha-fetoprotein
- Relaxin
- Distemper
- Parvovirus
- Rabies Virus
- Ehrlichia
- Toxoplasma Gondii
- Parvovirus/Coronavirus combined test
- Cystatin C (SDMA alternative)
- Anaplasma Phagocytophilum
- NT-proBNP/Troponin I combined test
- Distemper/Parvovirus/Infectious Hepatitis Virus combined test
- cCRP
- D-Dimer
- HbA1c
- NT-proBNP
- Troponin I
- HP (Helicobacter Pylori)
- Lyme Disease
- Coronavirus
- Giardia
- Rotavirus
- SDMA



LOW COST, PORTABLE CANINE IMMUNOASSAY TESTING

The InSight Canine-IA is a portable, handheld, easy to use Canine Immunoassay Analyser providing quick, accurate and reliable results.

FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Portable, with built-in rechargeable battery - ready to use anytime, anywhere
- Easy to use with touch screen technology
- Wi-Fi and Bluetooth connectivity - easy transfer of data
- Accurate and reliable
- Multi-language operation system
- Wi-Fi connectivity
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



Software Languages

English, Chinese, Spanish, French, Russian, Portuguese, Turkish, Italian, Polish, German

SPECIFICATIONS

Operating System:	Built-in Android test management software
Sample Type:	Serum and Plasma (EDTA)
Data Storage:	ROM Capacity: 8GB RAM Capacity: 1GB
Built-in Battery:	Removable rechargeable battery Standby Time: >72 hours Operating Time: >4 hours
Data Output:	USB data exporting Built-in thermal printer Data Transfer
Wireless Communication:	4G, Wi-Fi, Bluetooth
Image Processing:	Built-in 500W autofocus camera
Working Environment:	Temperature: +10°C to +30°C Humidity: ≤80% RH
Input Voltage:	DC 12V/5A
Dimensions:	240 (W) x 90 (D) x 65 (H) mm
Weight:	0.8kg

TESTS AVAILABLE

- | | |
|---|----------------------------|
| • cPL | • cCRP |
| • T4 | • D-Dimer |
| • TSH | • HbA1c |
| • Cortisol | • NT-proBNP |
| • Progesterone | • Troponin I |
| • Alpha-fetoprotein | • HP (Helicobacter Pylori) |
| • Relaxin | • Lyme Disease |
| • Distemper | • Coronavirus |
| • Parvovirus | • Giardia |
| • Rabies Virus | • Rotavirus |
| • Ehrlichia | • SDMA |
| • Toxoplasma Gondii | |
| • Parvovirus/Coronavirus combined test | |
| • Cystatin C (SDMA alternative) | |
| • Anaplasma Phagocytophilum | |
| • NT-proBNP/Troponin I combined test | |
| • Distemper/Parvovirus/Infectious Hepatitis Virus combined test | |



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.