

InSight[®] Crea-Vet

Veterinary Creatinine Meter



SIMPLE, ACCURATE
ASSESSMENT OF KIDNEY FUNCTION
IN 30 SECONDS.



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight Crea-Vet is a portable, handheld Veterinary Creatinine Meter that provides accurate, quantitative creatinine results in 30 seconds.

Creatinine measurements are used in the diagnosis, treatment and monitoring of renal diseases.

FEATURES

- Rapid assessment of kidney function
- Easy to use
- Portable
- Can be used in the field
- No calibration required
- Accurate, quantitative creatinine results in 30 seconds
- Small $\leq 2\mu\text{l}$ sample size
- Capillary or venous whole blood
- One-click strip ejection prevents cross contamination
- Stores 500 results and 20 QC results
- Bluetooth connectivity

SPECIFICATIONS

Unit of Measurement: $\mu\text{mol/L}$ or mg/dL

Sample Type: Fresh whole blood

Sample Volume: $\leq 2\mu\text{l}$

Detection Time: 30 seconds

Memory Capacity: 500 results and 20 QC results

Automatic Shutdown: After 2 minutes of inactivity

Power Supply: DC 3 V, 10 mA

Dimensions: 107 x 54 x 23 mm (± 5 mm)

Weight: 100 ± 15 g (including batteries)

Operating Temperature: $+15^{\circ}\text{C}$ to $+40^{\circ}\text{C}$

Operating Relative Humidity: 10% to 85% (non-condensing)

Atmospheric Pressure: 68 kPa to 106.6 kPa

Storage Temperature: -20°C to $+55^{\circ}\text{C}$

Storage Relative Humidity: 10% to 93% (non-condensing)



Simple, accurate assessment
of kidney function in 30 seconds



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.