

RapidVet[®] -H

Blood Typing and Crossmatch Kits - Canine and Feline

"Every Dog & Cat in need of a transfusion should be blood typed"



RAPID IN-PRACTICE
BLOOD GROUP TYPING
AND CROSSMATCH TESTS

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

RapidVet® -H

Blood Typing and Crossmatch Kits

RapidVet® - H Blood Typing Kits - Canine & Feline

In-practice tests for rapid determination of Canine DEA 1.1 and Feline A, B or AB blood types.

These fast in-house blood typing card kits for Dogs and Cats are essential for transfusion and breeding decisions. They are fast, easy to use, economical and can be stored at room temperature.


RapidVet-H Canine

- For Canine DEA 1.1
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



RapidVet-H Feline

- For Feline A, B & AB
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



"Every Dog & Cat in need of a transfusion should be blood typed" Dr. Urs Giger, Dr.med.vet., Dip ACVIM, ECVIM

RapidVet® - H Canine DAL

Canine blood types other than DEA 1 can also be important in transfusion decisions. RapidVet-H Canine DAL is intended for use to classify dogs as DAL positive or negative.

This information can be essential in donor selection, especially in the Dalmatian, Doberman, Shih Tzu and Lhasa Apso breeds, which reportedly have a higher percentage of DAL-negative blood types. It is also useful in determining the cause of an incompatible crossmatch between DEA 1 compatible dogs.

As many as 98% of dogs, including blood donors, carry the DAL antigen. This antigen elicits strong immunogenicity. Sensitisation through transfusion can occur in dogs that lack it.

DAL-negative blood types are more common in certain breeds and they are more at risk of delayed and acute haemolytic transfusion reactions if transfused with DAL positive blood. Doberman Pinschers are especially at risk since the prevalence of both DAL-negative blood and von Willebrand disease is relatively high.



- Enable more complete, in-practice, pre-transfusion work-up and donor selection
- Useful in determining possible cause of an incompatible crossmatch between DEA 1 compatible dogs
- Clear, easy to follow instructions
- Visible results within 1-2 minutes
- 2 or 5 card kits available

RapidVet® - H Crossmatch Kits Use alongside Canine & Feline Blood Typing Kits



In-practice crossmatch test kits for use in transfusion medicine on either Canine or Feline patients.

- Will detect all major Canine and Feline blood group incompatibilities
- Performed in-practice in 10 minutes
- Easy to use and results are clear
- Available in kits of 3 tests
- Minor Crossmatch Kits also available

**Essential for
good transfusion
practice**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**