

# Vortex Mixer

Available for hire or purchase for your clinical trial.



## Overview

This vortex mixer can be used for gentle mixing to vigorous agitation. Choose between 'touch' mode or continuous operation simply by pressing a button on the fascia.

Robust die-cast body avoids unnecessary movement during use. It has an ergonomic low profile design, making everyday vortexing comfortable for the user.

Variable speed: 200 - 2500 rpm

## Operating Feature

- Integral retort rod fixing allows vessels to be secured above the vortex action for long-term mixing
- Intermittent or continuous mode
- Stable low profile body
- Robust die-cast construction
- Variable speed control
- 'Touch' mode
- Can mix with a variety of other vessels when used with the accessory pack
- Accessory Pack (optional):
  - plastic cradle for standard microtitre plate
  - circular foam insert for beakers and flasks (up to 500ml)
  - rectangular foam insert (8 x 0.2ml, 8 x 0.5ml and 16 x 1.5ml microcentrifuge tubes)
  - solid rectangular foam insert (blank for custom holes)

## Specification

- Speed: 200 to 2500 rpm
- Orbit Diameter: 4.2mm
- Power Supply: 90 - 240V, 50/60 Hz, 20 W
- IP Rating: 31
- Weight: 3.2 kg
- Dimensions: 78h x 135w x 215d mm



**Please contact us for a quotation for the hire or purchase of this product**

Product description - Vortex Mixer

Product code - WD7VOR002

Product specifications can be subject to change without prior notice.

Version: V2 10/19

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company